

ScienTempTM

3565 S. Adrian Hwy
Adrian, MI 49221
Ph: 517-263-6020
www.scientemp.com

Operating and Installation for Super and Ultra Cold Countertop Freezers

For Model's: 45-01 86-01

Customer Service / Replacement Parts:

Email: tracy@scientemp.com
Ph: 517-263-6020 Ext 102

Technical Support:

Email: stuart@scientemp.com
Ph: 517-263-6020 Ext 101




TABLE OF CONTENTS

CONTACTING THE FACTORY.....	1
IMPORTANT BEFORE YOU START.....	1
PRECAUTIONS AND SAFETY.....	2
CHOOSING THE BEST LOCATION.....	3
START UP.....	4
GET TO KNOW YOUR FREEZER.....	4
CAL9500 DIGITAL TEMPERATURE CONTROLLER.....	4-5-6
LID SEAL.....	7
CHART RECORDER (IF EQUIPPED).....	7-8
MAINTENANCE.....	10-11
TROUBLESHOOTING.....	12-13
WARRANTIES.....	13
WIRING DIAGRAM.....	14

CONTACTING THE FACTORY

For reference and when contacting the factory, please have your freezer information ready:

Model Number:

Serial Number:

Date Received:

*The Model Number and Serial Number can be found on the data plate attached to the cabinet on the back upper left corner.

IMPORTANT...READ FIRST

- It's recommended to keep the freezer running at a temperature below freezing rather than turning the freezer off and on.
- Make sure your freezer is placed on a level surface and in a well ventilated room. Provide at least 3" of space around the freezer for ventilation. Avoid placing the freezer in locations exposed to direct sunlight or any heat source.
- Alarm On/Off switch should remain in the "Off" position until freezer reaches desired set point temperature. When the temperature has been reached, turn the alarm switch to the "On" position.
- Freezers are pre-programmed for proper operation.
**Depending on your specifications, you may need to only adjust the set point.
- Freezer should not be loaded with product until the cabinet has operated for 24 hours at the desired setpoint temperature.
- Caution: When skin is exposed to low temperature conditions, frostbite is possible.
**Please wear proper gear to avoid freezing of exposed skin.
- The outer surface should be kept clean and free of any wrapping, covers or objects that will limit the transfer of the heat from the freezer shell.
- If loading large amounts of product into the freezer at one time, be aware that the inside temperature will increase and may trigger the alarm. (You may turn the alarm off until the unit has fully recovered back to set point temperature)
- This freezer is not a *rapid-freeze* unit.

PRECAUTIONS AND SAFETY

In this manual, the words WARNING and CAUTION mean the following:

- *  **WARNING** A potentially hazardous situation, which if not avoided, could result in serious injury or death.
- *  **CAUTION** A potentially hazardous situation, which if not avoided, may result in minor or moderate injury or damage to equipment.

Before installing, using or maintaining this product, please be sure to read this manual carefully. Failure to follow these instructions may cause this product to malfunction, which could result in injury or damage.

Below are important precautions that apply to this product:

- Use this product only in the manner described in this manual. Before using it, verify that this product is suitable for its intended use.
- Do not modify the cabinet by drilling new holes.
- The cabinet must be properly grounded in accordance with the national and local electrical codes and ordinances.
- Disconnect the unit from all power sources before doing any maintenance or cleaning of the unit.
- If you are storing high value product in freezer, it is recommended that an independent monitoring system be used to monitor and alert you of any issues that may interfere with the freezer's operation.

CHOOSING THE BEST LOCATION

MOVING YOUR FREEZER:

The refrigeration system is designed to operate with the cabinet located on a flat surface.

****Do not tilt the cabinet more than 30° to any side.****

If the cabinet must be tilted on an angle for handling or moving purposes, allow it to sit in an upright position for at least two (2) hours prior to plugging in and starting the freezer.

THE LOCATION:

Wiring and Installation:

Your freezer is designed for indoor use only. Please confirm your power supply is adequate and delivers the correct voltage requirements for your new freezer. Operating the freezer with insufficient power can result in severe damage to the compressor and electrical components.

~Voltage requirements for your freezer can be found on the serial tag on back upper left side of your cabinet~

115V cabinets are equipped with a Nema 5-15P three-prong plug for your protection against shock hazards.

Choose a location that is close to a dedicated electrical line with a proper grounded (three-prong) electrical outlet.

****Note:** 230V freezers do not come with a plug, customer is responsible for installing one according to electrical configuration.

⚠️WARNING ***Do not use an extension cord or surge protector.**

⚠️WARNING ***Do not use an adapter of any sorts.**

⚠️WARNING ***Do not remove or modify the ground pin on the electrical plug.**

⚠️WARNING ***To avoid freezer being shut off....Do not plug freezer into outlet controlled by a wall switch.**

⚠️WARNING ***When unplugging the unit from the wall outlet, never pull on the cord. Always grip the plug safely and pull straight out from the outlet.**

Air Circulation:

The cabinet should be situated in a well ventilated room to allow proper air circulation. Provide at least 3" of space around the entire freezer for ventilation.

***NOTE:** The cabinet must be installed on a sturdy, solid, level floor, table or countertop resting firmly on all four mounting points*

Ambient Temperature:

Do not place the freezer in locations exposed to direct sunlight, heat registers or any other heat sources.

Drafts from fans, air conditioning or open doors can also affect the cabinet and product temperature.

The ambient temperature range for the location must be 59° to 90°F (15° to 32°C) for optimal freezer performance.

***When opening the door, avoid leaving it open for extended periods of time. Chamber air will escape rapidly.**

Room air, which is higher in humidity, will replace chamber air and may cause frost to develop in the chamber more rapidly

START UP AND GET TO KNOW YOUR FREEZER

START UP:

1) Plug in or connect freezer power source

Interior Compartment Temperature:

The temperature is controlled by a CAL 9500 digital temperature controller. When you first power up the freezer, the high stage compressor and fan will start and operate for 3 to 6 minutes before the low stage compressor freezer, starts. Pull down to minimum temperature will take from 3 to 6 hours, depending on ambient temperature and size of cabinet. The high stage compressor runs 100% of the time. The low stage compressor will cycle on and off to maintain the temperature of the controller setting. Set the control to the desired set point.

2) Set your set point temperature:

To view the set point temperature, press the * key. Press and hold the * key and use the up or down arrow to increase or decrease the set point temperature.

Freezers are shipped from the factory with the indicating digital control programmed for proper operation.

INSTRUMENT PANEL FEATURES

! This page can be photocopied and used as a visual aid and bookmark when working in other parts of the manual.

Green LED:

Setpoint 1 output indicator

Green Display:

Process variable or
Function/Option

Upper Red LED:

Setpoint 2 output
indicator

Lower Red LED:

Setpoint 3 output
indicator



Right Red LED:

Program Holdback
indicator

Orange Display:

Setpoint value or
program selection

ADJUSTMENTS

To enter or exit **program mode**:

To scroll through **functions**:

To change **levels** or **options**:

To view setpoint units:

To increase setpoint:

To decrease setpoint:

To reset latched alarm or tune fail:

To run or Hold a program:

Press ▲ ▼ together for 3 seconds

Press ▲ or ▼

Press * ▲ together or * ▼ together

Press *

Press * ▲ together

Press * ▼ together

Press ▲ ▼ together briefly

Press * ▼ together for 3 seconds

Notes:

If in difficulty by becoming "lost" in program mode, press ▲ and ▼ together for 3 seconds to return to display mode, check the INSTRUMENT ADJUSTMENTS above and try again.

When in program mode, after 60 seconds of key inactivity the display will revert to either *inPt : nonE* or, if the initial configuration has been completed, the measured value. Any settings already completed will be retained. During Program Configuration it is recommended that this feature is inhibited. Select **ProG StAY** in Level 4.

CAL9500 Configuration Data

****Note: Your CAL9500 Controller is preset to the following settings. No need for programming. This page is for informational purpose only.**

Level 3

SP1.d – rly

SP2.d – rly

Burn – uP.SC

rEV.d – 1d.2r

rEV.L – 1i2n

SPAn – 0.0

*Zero – 3 *Note: only change this parameter if calibrating unit, value may vary

ChEk – off

rEAd - Var

tECh – CtA

Ver – 953

rSET – none

Level 2

SP1.P – 0

hAnd – off

PL.1 – 100

PL.2 – 100

SP2.A – Dv.hi

SP2.b – none

di SP – 0.1

hi .SC – 100

Lo. SC – -100

in Pt – tc t

unit – C or F

Level 1

tune – off

band – 0.3

int.t – off

dEr.t – off

dac – 0.5

cyc.t – on.of

ofSt – 0.0

SP.Lk – off

SPrr – 0

SPrn – off

SoAk – --

Set.2 – 8 (alarm deviation)

bnd.2 – 0.1

Cyc.2 – on.of

English

English

Range of Adjustment shown under description.

If applicable, factory settings shown in **bold**.

For a full description of menu functions refer to pages 7 to 10, and 16/17 in Programmer section.

Note: The letter K appears in the instrument display as the character K

Lid Seal:

The magnetic door gasket needs to seal firmly on all four sides. Remove any excessive frost build up and keep the gasket clean. Failure to remove the frost regularly could result in damage to the door and/or gasket along with not allowing it to properly close.

Alarm Signal System:


On start up: The Alarm should be kept in the “OFF” position until the cabinet set point temperature has been reached. When the temperature has been reached, switch to the “ON” position.

The Alarm is equipped with a time delay relay to prevent the alarm from going off every time the door is opened. This relay can be adjusted from 1 to 10 minutes. This delay is factory set for 10 minutes.

To adjust the time delay relay:

- a. Remove screw located on right side, middle of the instrument panel with the temperature controller. Swing instrument door open. Time delay has a dial and is located on lower right corner of electric box.
- b. If you have a chart recorder: Open the glass door to get to the black instrument panel with the chart recorder. Remove screw located on right side, middle of panel. Instrument panel will hinge open. Time delay has a dial and is located on lower right corner of electric box.

Chart Recorder (If equipped) :



4" Chart Recorders are an optional item for countertop freezer models. Recorders are used for recording the temperature inside the freezer.

NOTE: On start up you must replace the round paper chart, connect the 9V battery strap to battery and remove the black rubber cap on ink pen. Recorder will start as soon as freezer is connected to power. But may not start to record temperature until freezer reaches the temperature range that the recorder is set up for.

****While you have the arm moved off the paper for changing the chart, remove the rubber cap.**

Getting started with your recorder:



- 1) Open cover to access the recorder.
- 2) Make sure recorder has new 4" recorder paper ready for use. If you need to replace the paper, please follow the instructions under "Changing the Chart Paper" section below.

Changing the Chart Paper:

- 1) Locate "Change Chart" key pad on top left corner on recorder.
 - 2) Press and hold the Change Chart button. (#3) for approximately 1 second until the pen begins to move to the left and then release the button. When the pen has moved completely off of the chart, unscrew (counter-clockwise) the chart "hub" knob located in the center of the paper. Remove the old chart. When installing the new chart, make sure you position the new one so that the correct time line coincides with the time line groove on the chart plate. See Diagram on next page to help with location of the time line groove.
 - 3) After you have the new paper in place, re-attach the "hub" knob and make sure its secured in place by hand. ****Do not over tighten. Over tightening can damage the gears in the chart motor gearbox assembly.**
 - 4) Press #3 button again for about 1 second until pen starts to move back onto the chart and then release the button.
- NOTE:** If red ink pen is not marking on paper, you may need to replace it.

Changing the Red Fiber Tipped Pen:

- 1) For ease of replacment, it is suggested that the two screws that hold the metal pen arm be loosened and the pen cartridge and metal pen arm be removed as an assembly. ****See diagram on prior page for location of pen arm screws.**
 - 2) Unsnap the plastic "U" clip tab of the pen cartridge from the metal pen arm, remove and discard the old pen cartridge.
- Replace the new cartridge by opening the hinge and snapping it securely around the metal pen arm.



Pen Arm Calibration :

To check and/or adjust the recording pen calibration to the outer most temperature graduation of the chart, press and hold the "change chart" button (#3) until the pen begins to move off of the chart. Once the pen has moved off of the chart, again press and hold the "Change Chart" button (#3) until the pen begins to move back onto the chart. The pen should briefly stop at the outer most temperature graduation of the chart before continuing onto the chart to begin recording. If the pen does not stop exactly at this location on the chart, it can be adjusted or "calibrated" by using the left (#1) or right (#2) arrow buttons.

When the pen moves back onto the chart and briefly stops, you will have approximately 5 seconds in which to adjust the pen's position using the left and right arrow buttons.

Temperature Recorder Calibration Check :

NOTE: The recorder has been accurately calibrated at the factory. Before making any adjustments, the freezer with recorder should be in service and running for at least 24 hours.

If any adjustments are needed, perform the following steps:

- 1) Place a Certified Test Thermometer in a solution bottle inside of the freezer alongside or nearest to the recorder's sensor probe.
- 2) Once the freezer and solution level out, compare the position of the pen on the recorder to the test thermometer's reading
- 3) If an adjustment is needed, using the left (#1) and right (#2) arrow, Push button on the recorder to calibrate (or move) the pen's position on the chart to correspond to the temperature of the solution.

NOTE: The arrow buttons must be pressed for approximately five (5) seconds before the pen will begin to move.

MAINTENANCE

* **CAUTION** Maintenance or repair of your freezer should only be performed by a trained professional.

*The refrigeration system is hermetically sealed and the fan motor has lifetime lubrication, they requires no maintenance.

~It's recommended to keep the freezer operating at a temperature below freezing rather than turning the freezer off and on.

CLEANING THE CONDENSER:

The finned condenser in the refrigeration compartment of the freezer, should be cleaned of dust and debris about once a year (or more often if the freezer is located where excessive dust and dirt may occur)

To clean the condenser:

- Remove power to freezer.
- Use forced air and a vacuum to clean dirt/debris from the condenser and refrigeration compartment.

How to clean your condenser on your 45-01/86-01

Depending on your model you will need to first remove a screw.




If your model has a Chart Recorder

#1
Remove this screw



If your model does not have a Chart Recorder

Open up to the compartment behind and you will see two more screws



#2
Remove these screws

Next, gently separate and slide the right grill panel off. (DO NOT remove the left side)



#3



#4

Slide this part completely off to the right.

#5 - After the grill is removed you will see the black condenser. Being careful not to bend the steel fins, vacuum vertically up and down.

MAINTENANCE

DEFROSTING THE FREEZER:

Defrosting: Accumulation of frost may be removed in a couple of ways.

1.) While the freezer is running the frost may be removed with a plastic or wooden scraper or spatula, do not use metal or sharp objects such as an ice pick. This will prevent serious damage to the finish of the inner-liner of the freezer. The frost and ice may be collected in a pan, towel or any device to prevent it from falling on the floor or into the freezer onto the material stored inside the freezer.

2.) In the case where the freezer is not loaded, frost and ice may be removed by disconnecting the power the power until the frost and ice thaw. After defrosting, dry the freezer compartment interior and plug back in. **To help with the defrosting procedure, open lid to allow warm air inside the chamber.

MAINTENANCE FOR THE LID GASKET:

Occasionally check and make sure there is no frost or ice build up on or around the gasket. Keeping the gasket clean and free of ice and frost will help prevent the gasket from getting torn or punctured. Make sure the gasket is completely connected on all sides and laying flush to the freezer surface.

CLEANING THE EXTERIOR AND INTERIOR:

As needed: Use only a soft cloth and mild cleaners. Do not use harsh solvents or abrasives to clean the interior or exterior of the freezer.

CHANGING BATTERIES:

It's recommended to check your accessory batteries regularly and replace if there are signs of corrosion or leakage. Yearly replacement is recommended to ensure your equipment is functioning properly and ready in an event of power failure. Batteries can be found on your freezer in two areas.

1) Chart Recorder - (If equipped) (1) 9V Battery - Found on top right corner of the Recorder.

TROUBLESHOOTING GUIDE

PROBLEM	POSSIBLE CAUSE	SOLUTION
Freezer does not power up / no digital display	1. Freezer not plugged into power supply	1. Plug freezer into dedicated power source (wall outlet)
	2. Power cord not firmly seated in power outlet.	2. Push plug firmly into outlet
	3. No power at outlet	3. Move freezer to different outlet or try plugging something else into outlet to test the outlet.
	4. Fuse Blown or Circuit Breaker tripped	4. Check breaker/ fuse box and make sure that it's not tripped or have blown fuse.
	5. Make sure you have correct voltage at outlet	5. Contact service technician
Freezer is not cooling / Digital Display is on	1. Compressor not running	1. Possibly bad compressor or compressor start components. (Contact service technician)
	2. Possible leak or no refrigerant	2. Contact service technician
	3. Check voltage	3. Make sure there is proper voltage at outlet. (Contact service technician)
	4. Check for proper controller settings	See "CAL9500 Configuration Data" for proper settings
Temperature inside freezer does not reach the desired setpoint temperature / recovers slowly	1. Condenser coils are plugged or dirty	1. Contact service technician to properly clean coil
	2. Freezer Gasket is dirty or damaged	2. Clean or replace the bad gasket
	3. Ambient temperature is too warm	3. Cool location down and keep away from heat source
	4. Large amounts of warm product was added or too much product	4. Allow Freezer to recover after loading with warm product.
	5. Door may be opened to frequently	5. Warm air entering freezer will cause the freezer to warm up.
	6. Not enough air circulating around freezer	6. Make sure there is at least 3" of space around the outside of the freezer.
	7. Make sure condenser fan motor is operating/spinning in refrigeration compartment.	7. If fan motor is not spinning, determine the reason and fix.
Freezer too cold	1. Controller is set too cold	1. Change temperature
	2. Power relay is faulty	2. Contact service technician
External surface of freezer is warm/hot to touch	1. Warmer than room temperature	1. This is normal at start up. Freezer compressor(s) transfers heat from inside the unit to the outside, through the walls.
Excessive frost build up	1. Hot and humid air	1. This is normal.
	2. Door opened frequently or long periods of time	2. Keeping door closed helps with frost build up.
	3. Door may not be sealing correctly	3. Check gasket area for frost build up. Remove any frost and keep gasket clean. *Also check for damage to the gasket.
		~ Make sure the lid is flat and seal is seated correctly
Door will not close	1. Excessive frost build up	1. Using a wood scraper carefully remove excessive frost

TROUBLESHOOTING GUIDE

Compressor will not start - no hum	1. Power cord not plugged in at power source.	1. Plug freezer into power outlet
	2. Wiring improper / loose / disconnected	2. Check wiring
	3. No power to outlet	3. Make sure power outlet has power being supplied to it.
		~Check breaker/ fuse box and make sure that it's not tripped or have blown fuse.
Compressor will not start - hums but trips on overload protector	1. Low voltage to unit	1. Determine the reason and correct.
	2. Starting capacitor defective.	2. Determine the reason and replace.
	3. Compressor motor has a winding open or shorted	3. Determine the reason and correct, replace if necessary.
	4. Extension Cords	4. Do not use extension cord. Must be plugged directly into wall outlet using power cord attached to freezer.

NOTES:

[illegible]

US / CANADA WARRANTY



Scientemp

3565 S. Adrian Hwy ♦ Adrian, MI 49221
800-968-2653 ♦ 517-263-6020 ♦ Fax 517-263-5492

Standard Warranty

Applies Only Within United States and Canada

SCIENTEMP WARRANTS TO THE ORIGINAL PURCHASER

FIRST 18 MONTHS – The cabinet and all of its parts, shall be free of defects in material and workmanship under normal use and service for a period of 18 months from the date the unit has been shipped from our facility. Scientemp's sole obligation under this warranty shall be limited, at its option, to either repairing or replacing any part of the cabinet determined by an authorized service agent to be defective. Scientemp reserves the right to repair the freezer at our facilities.

THE SCIENTEMP STANDARD WARRANTY DOES NOT COVER

TRANSPORTATION COSTS – Scientemp shall not be responsible for transportation or incidental costs incurred in connection with the repair or replacement of a cabinet or any of its parts.

ABUSE, MISUSE, ACCIDENTS – Scientemp shall not be responsible for parts or assemblies which upon inspection are determined by an authorized Scientemp Service Agent to have been subjected to misuse, neglect, alteration, accident, abuse, damage during transit or delivery, or by fire or flood.

CONSEQUENTIAL DAMAGES – IN NO EVENT SHALL SCIENTEMP CORP BE LIABLE FOR SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES, INCLUDING WITHOUT LIMITATION, PRODUCT LOSS OR PRODUCT SPOILAGE CLAIMS, NOR FOR ANY DELAY IN THE PERFORMANCE OF THIS WARRANTY DUE TO CAUSES BEYOND ITS CONTROL.

GENERAL

The standard warranty and any service contract related to the STANDARD WARRANTY shall apply only to the products sold and used within the boundaries of the United States and Canada.

Users may file warranty claims either directly with Scientemp Corp, 3565 S. Adrian Hwy, Adrian, MI 49221, or with the seller from whom the cabinet was purchased. All claims must be supported by information concerning the alleged defect and specifically identified by the Serial Number of the cabinet.

THERE ARE NO OTHER WARRANTIES EXPRESS, IMPLIED, OR STATUTORY, EXCEPT THIS WARRANTY, WHICH IS IN LIEU OF ALL OTHER WARRANTIES INCLUDING TO THE EXTENT PERMITTED BY LAW, ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

DO NOT DRILL HOLES IN CABINET

Refrigeration tubing and wiring is routed through the cabinet walls. Leaks, wet insulation or electrical problems caused by drilling holes are not covered by warranty.

INTERNATIONAL WARRANTY



Scientemp

3565 S. Adrian Hwy ♦ Adrian, MI 49221
800-968-2653 ♦ 517-263-6020 ♦ Fax 517-263-5492

International Warranty

SCIENTEMP WARRANTS TO THE ORIGINAL PURCHASER

FIRST 18 MONTHS – The cabinet and all of its parts shall be free of defects in material and workmanship under normal use and service for a period of 18 months from the date the unit has been shipped from our facility. Under this warranty, Scientemp will supply replacement parts, F.O.B. Adrian, Michigan, U.S.A., for any part proven defective in material and workmanship. Scientemp is not responsible for supplying the labor required to replace parts and make necessary repairs.

THE SCIENTEMP STANDARD WARRANTY DOES NOT COVER

TRANSPORTATION COSTS – Scientemp shall not be responsible for transportation or incidental costs incurred in connection with the repair or replacement of a cabinet or any of its parts.

ABUSE, MISUSE, ACCIDENTS – Scientemp shall not be responsible for parts or assemblies which upon inspection are determined by an authorized Scientemp Service Agent to have been subjected to misuse, neglect, alteration, accident, abuse, damage during transit or delivery, or by fire or flood.

CONSEQUENTIAL DAMAGES – IN NO EVENT SHALL SCIENTEMP CORP BE LIABLE FOR SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES, INCLUDING WITHOUT LIMITATION, PRODUCT LOSS OR PRODUCT SPOILAGE CLAIMS, NOR FOR ANY DELAY IN THE PERFORMANCE OF THIS WARRANTY DUE TO CAUSES BEYOND ITS CONTROL.

GENERAL

The standard warranty and any service contract related to the STANDARD WARRANTY shall apply only to the products sold and used outside the boundaries of the Continental United States.


Users may file warranty claims either directly with Scientemp Corp, 3565 S. Adrian Hwy, Adrian, MI 49221, or with the seller from whom the cabinet was purchased. All claims must be supported by information concerning the alleged defect and specifically identified by the Serial Number of the cabinet.

THERE ARE NO OTHER WARRANTIES EXPRESS, IMPLIED, OR STATUTORY, EXCEPT THIS WARRANTY, WHICH IS IN LIEU OF ALL OTHER WARRANTIES INCLUDING TO THE EXTENT PERMITTED BY LAW, ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

DO NOT DRILL HOLES IN CABINET

Refrigeration tubing and wiring is routed through the cabinet walls. Leaks, wet insulation or electrical problems caused by drilling holes are not covered by warranty.

45-01 WIRING DIAGRAM - CAL 9500



86-01 WIRING DIAGRAM - CAL 9500

